



# FISH™ CONTROL CLOSE™

Maintain access and control while you perform vascular closure.

ControlClose™ enables physicians to deliver healing to the artery following a catheterization procedure. It is able to be deployed regardless of vessel size, disease state or stick location. ControlClose™ also provides closure over-the-wire.

## Biologic Comfort

ControlClose™ and CombiClose™ heal the femoral artery by using a proven, comfortable, biologic material called SIS® (small intestinal submucosa).

\* Learn more about SIS® at [www.cookbiotech.com](http://www.cookbiotech.com)



# FISH™ COMBI CLOSE™

Perform your procedure and closure, over-the-wire, with the same naturally versatile device.

CombiClose™ can be used as a combined introducer sheath and closure device. It also provides patient comfort, ability to close any vessel, and over-the-wire deployment.

# CONTROL CLOSE™



Product Code	French Size	Hub Color
CS-6F-5-S	6 Fr	Green
CS-7F-5-S	7 Fr	Orange

Each device includes one FISH™ ControlClose™ Sheath w/ Closure Patch and Stabilizer, Dilator and Guide Wire. FISH™ ControlClose™ devices are individually packaged and sold in boxes of (5) devices.

# COMBI CLOSE™



Product Code	French Size	Hub Color
16-05-05-S	5 Fr	Gray
16-06-05-S	6 Fr	Green
16-07-05-S	7 Fr	Orange
16-08-05-S	8 Fr	Blue



Each device includes one FISH™ CombiClose™ Sheath w/ Closure Patch, Dilator, Guide Wire, Pre-Dilator and an Obturator. FISH™ CombiClose™ devices are individually packaged and sold in boxes of (5) devices.

## CLINICAL EVIDENCE

**Complete healing and no loss of future access at puncture site.**

The 6-week post procedure image (A) demonstrates complete healing and no loss of future access at the original puncture site. This healing provides increased patient comfort.

Both devices are capable of closing arteries as small as 3mm. Image (B) depicts an artery that was successfully closed with a FISH™ device.



[www.morrisinnovative.com](http://www.morrisinnovative.com)

Morris Innovative, Inc.  
907 West Second Street  
Bloomington, IN 47403  
812 355 0450  
509 356 7987 fax

Affiliate / Distributor Information

